


Table Top pH Meter

Low
Cost

WITH
A T C



Model : Equinox-pH 101A

EquinoxTM

Economy Bench Top pH Meter

1 Features

Bright Red LED	: Better Visibility even in 'Hard - to See' areas.
Auto Temp Compensation(ATC)	: Built-In
Accuracy	: Good accuracy levels to match the Industry Standard.
Calibration	: 2 Point calibration for better stability and accurate readings.
Electrode Holder	: Provided, for added safety to the electrode.
Warranty	: Unprecedented 3 year warranty on the Instrument.

2 Specifications

Range : 0.00 to 14.00 pH ±1999 mV 0 to 100 °C	Temperature : 0 to 100 °C Auto Compensation (using RTD Sensor)
Resolution : 0.01 pH 1 mV 0.1 °C	Analog output : 4-20 mA for 0.00 to 14.00 pH
Accuracy : ± 0.05 pH (without Electrode) ± 2mV ± 0.5 °C	Calibration : 2 points
Input Impedance : 10 ¹² Ohm	Connection : BNC
	Power Supply : 230 VAC ±10%, 50Hz
	Size : 220 (W) × 90 (H) ×175 (D)mm
	Weight : Approx 2 Kg.

3 Ordering Information

Order No.	Item
Equinox-pH101AC	Complete Kit consisting of pH Meter , ATC Sensor, Glass pH Electrode, Electrode Stand with clip, One Bottle (200 ml) each of pH Buffers 7.0, 4.0, 10.0
Accessories :	
Equinox-pH101A	Consists only pH meter, Electrode stand with clip
Equinox-pH101A / 01R	ATC sensor with 1 meter cable
Equinox-pH101A / 01E	Glass pH electrode
Equinox-pH101A / 01S	Electrode stand with clip
Equinox-pH101A / 01D	Low Ionic strength pH Electrode (DM water application)
Equinox-pH101A / 01EPX	Epoxy Body pH Electrode
Equinox-pH101A / 01B	Ph 7.0 buffer Solution (200 ml)
Equinox-pH101A / 02B	Ph4.0 Buffer Solution (200 ml)
Equinox-pH101A / 03B	pH 10.0 Buffer Solution (200 ml)

Authorised Dealer :



Valli Aqua and Process Instruments

1/231 A, Pananchalai Street, Adayalampattu,
(Near Raja Rajeshwari Engg. College)
Chennai - 600 095, Tamil Nadu.
Mobile : +91-9840752690
E-mail : enquiry@vapinstruments.in, kvgiri_59@yahoo.in
Website : www.vapinstruments.in

Note : We reserve the right to make changes, improvements and modifications to products shown.



EQUINOX