

TABLE TOP TDS METER

Model : Equinox - TDS 101

WITH
ATC



EquinoxTM

Economy Bench Top TDS Meter

An ergonomically designed lab model TDS meter housed in anti- corrosive enclosure

1 Features

Bright Red LED	:	Better Visibility even in 'Hard - to See' areas.
Auto Temp Compensation(ATC)	:	Built-In
Accuracy	:	Good accuracy levels to match the Industry Standard.
Calibration	:	1 Point calibration for better stability and accurate readings.
Electrode Holder	:	Provided, for added safety to the electrode.

2 Specifications

Measurement Range	: 0.0 to 199.9 ppm 0 to 1999 ppm 0.0 to 19.99 ppt 0 to 199.9 ppt	Temp. Compensation	: Auto 0 – 60°C
Accuracy	: 1% of Fs. (Excluding Rbbe Error)	Range Selection	: Manually Using R otor switch
Resolution	: 0.1 ppm 1 ppm 0.01 ppt 0.1 ppt	Power Supply	: 230 VAC 1 50 Hz
		Dimensions	: 230 (W)×180(D)×90(H)mm

3 Ordering Information

Order No.	Item
Equinox-TDS 101C	Complete Kit consisting of TDS Meter Conductivity Electrode of Glass body and of cell constant = 1, Electrode Stand with Clip, One Bottle (200 ml) TDS Buffer Standard of 1000 ppm
Accessories :	
Equinox-TDS 101	Consists only TDS Meter with electrode
Equinox- TDS101/01E	TDSElectrode of Glass Body and K=1
Equinox- TDS101/02E	TDSElectrode of Glass Body and K=0.1 (for low TDS Measurement)
Equinox- TDS101/01S	Electrode Stand with Clip
Equinox- TDS101/01B	TDSBuffer Standard of 1000 ppm
Equinox- TDS101/02B	TDSBuffer Standard of 7000 ppm

Authorised Dealer :



Valli Aqua and Process Instruments

1/231 A, Pananchalai Street, Adayalampattu,
(Near Raja Rajeshwari Engg. College)
Chennai - 600 095, Tamil Nadu.
Mobile : +91-9840752690
E-mail : enquiry@vapinstruments.in, kvigiri_59@yahoo.in
Website : www.vapinstruments.in

Note : We reserve the right to make changes, improvements and modifications to products shown.

TM
EQUINOX